

Product datasheet for **TA385930**

RTN4RL2 Rabbit Monoclonal Antibody [Clone ID: P14]

Product data:

Product Type:	Primary Antibodies
Clone Name:	P14
Applications:	ELISA, FC, IP, WB
Reactivity:	Mouse, Rat, Human
Host:	Rabbit
Isotype:	IgG, kappa
Clonality:	Monoclonal
Immunogen:	(His)6-tagged fragment of rat NgR2 (amino acids 279-420).
Specificity:	The antibody binds to human NgR2 with a Kd of 2.3 nM and to rat NgR2 with a Kd of 0.48 nM (bivalent binding). The antibody binds specifically to residues 351-399 of NgR2, a member of the Nogo-66 receptor family which inhibits axonal growth and sporulating in the CNS.
Formulation:	PBS with 0.02% Proclin 300.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid freeze and thaw cycles.
Stability:	3 years from dispatch.
Gene Name:	reticulon 4 receptor-like 2
Database Link:	Entrez Gene 349667 Human Q86UN3
Synonyms:	NgR2; NGRH1; NGRL3



[View online »](#)