

## Product datasheet for **TA385955**

### SIRPG Mouse Monoclonal Antibody [Clone ID: OX117]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	OX117
Applications:	ELISA, FC, IP
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1, kappa
Clonality:	Monoclonal
Immunogen:	Recombinant human SIRPy:CD4 fusion protein.
Specificity:	This antibody binds to human SIRP-beta and gamma.  This antibody binds to human Signal-regulatory protein gamma, an immunoglobulin-like cell surface receptor binding CD47. It mediates cell-cell adhesion and enhances T-cell activation.
Formulation:	PBS with 0.02% Proclin 300.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid freeze and thaw cycles.
Stability:	3 years from dispatch.
Gene Name:	signal regulatory protein gamma
Database Link:	<a href="#">Entrez Gene 55423 Human Q9P1W8</a>
Synonyms:	bA77C3.1; CD172g; FLJ42230; OTTHUMP00000185041; SIRP-B2; SIRP-beta-2; SIRPB2; SIRPgamma



[View online »](#)