

Product datasheet for **TA385966**

TNF alpha (TNF) Rabbit Monoclonal Antibody [Clone ID: CDP 571 (Humicade)]

Product data:

Product Type:	Primary Antibodies
Clone Name:	CDP 571 (Humicade)
Applications:	ELISA
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG, kappa
Clonality:	Monoclonal
Immunogen:	Human TNF-alpha
Specificity:	The antibody binds to human TNF-alpha with a Kd of 0.087 nM.

The antibody binds specifically to TNF-alpha, a pro-inflammatory cytokine which stimulates release of other cytokines e.g. IL-1, IL-6 and GM-CSF as well as secretion of ICAM-1 and VCAM-1 from endothelial cells. It stimulates acute phase responses and acts as a chemoattractant for NK cells and monocytes/macrophages. The antibody reached Phase II clinical trials with Celltech, UCB.

Formulation:	PBS with 0.02% Proclin 300.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid freeze and thaw cycles.
Stability:	3 years from dispatch.
Gene Name:	tumor necrosis factor
Database Link:	Entrez Gene 7124 Human P01375
Synonyms:	Cachectin; DIF; TNF-a; TNF-alpha; TNFA; TNFSF2
Note:	NOT FOR THERAPEUTIC USE - This is a research-grade biosimilar. This chimeric rabbit antibody was made using the variable domain sequences of the original Human IgG4 format, for improved compatibility with existing reagents, assays and techniques.



[View online »](#)