

## Product datasheet for **TA385975**

### dsRNA Mouse Monoclonal Antibody [Clone ID: 1D3]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	1D3
Applications:	ELISA, IHC, WB
Reactivity:	Virus
Host:	Mouse
Isotype:	IgG1, kappa
Clonality:	Monoclonal
Immunogen:	RDV-RNA-methylated bovine serum albumin complex.
Specificity:	

The antibody binds specifically to ds-RNA, which is part of the genetic material of some viruses, such as rotaviruses and the bluetongue virus. The antibody can be used as a diagnostic tool, to determine the nature of an unknown pathogen, and could potentially be used in gene therapy.

Formulation:	PBS with 0.02% Proclin 300.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid freeze and thaw cycles.
Stability:	3 years from dispatch.
Note:	This reformatted mouse antibody was made using the variable domain sequences of the original Mouse format, for improved compatibility with existing reagents, assays and techniques.



[View online »](#)