

Product datasheet for **TA385988**

LeY Mouse Monoclonal Antibody [Clone ID: H18A]

Product data:

Product Type:	Primary Antibodies
Clone Name:	H18A
Applications:	ELISA, IHC
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1, kappa
Clonality:	Monoclonal
Immunogen:	Lewis Y.
Specificity:	

The antibody binds specifically to LeY tetrasaccharide, a carbohydrate antigen expressed at the cell surface of selective carcinomas, including those of the colon, breast, ovary and lung. The monoclonal antibody may be of clinical use for the in vivo diagnosis and therapy of certain cancers.

Formulation:	PBS with 0.02% Proclin 300.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid freeze and thaw cycles.
Stability:	3 years from dispatch.



[View online »](#)