

## Product datasheet for **TA385990**

### IgE Mouse Monoclonal Antibody [Clone ID: TES-C21]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	TES-C21
Applications:	BI, ELISA
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1, kappa
Clonality:	Monoclonal
Immunogen:	Polyclonal human IgE.
Specificity:	

The antibody binds specifically to IgE, an immunoglobulin isotype secreted by upon stimulation of B cells by interleukins IL-4 and IL-13. IgE provides immunity to parasites such as helminths and *Plasmodium falciparum*. The immunoglobulin also plays a central role in the development of Type I hypersensitivity reactions and allergic diseases such as asthma, sinusitis and food allergies. These allergic responses are triggered by the binding of IgE to the high affinity receptor FcεRI, expressed on the surface of basophils and mast cells. Binding of IgE causes cross-linking of the receptor, inducing degranulation and release of inflammatory mediators, particularly of histamine. The monoclonal antibody TES-C21 binds to secreted IgE and to IgE expressed on the surface of IgE-producing B cells, thereby inhibiting IgE binding to FcεRI and the release of histamine from mast cells. The monoclonal antibody can have diagnostic applications, such as for identifying and quantifying IgE-expressing B cells and IgE levels in serum samples, and therapeutic applications in the treatment of allergies. Furthermore, the antibody can be conjugated to cytotoxic drugs and deliver them specifically to IgE-producing B cells.

Formulation:	PBS with 0.02% Proclin 300.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid freeze and thaw cycles.



[View online »](#)

**Stability:** 3 years from dispatch.

**Note:** This reformatted mouse antibody was made using the variable domain sequences of the original Mouse IgG format, for improved compatibility with existing reagents, assays and techniques.