

## Product datasheet for **TA385999**

### IL7R alpha (IL7R) Rabbit Monoclonal Antibody [Clone ID: YIL7R323]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	YIL7R323
Applications:	ELISA, FC
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG, kappa
Clonality:	Monoclonal
Specificity:	This antibody binds to an epitope of human Interleukin 7 receptor distinct to that bound by Ab00531.
Formulation:	PBS with 0.02% Proclin 300.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid freeze and thaw cycles.
Stability:	3 years from dispatch.
Gene Name:	interleukin 7 receptor
Database Link:	<a href="#">Entrez Gene 3575 Human P16871</a>
Synonyms:	CD127; CDW127; IL-7R-alpha; IL7RA; ILRA
Note:	This chimeric rabbit antibody was made using the variable domain sequences of the original Rat IgG format, for improved compatibility with existing reagents, assays and techniques.



[View online »](#)