

## Product datasheet for **TA386003**

### **Cd8b1 Rat Monoclonal Antibody [Clone ID: YTS 156.7.7]**

#### **Product data:**

<b>Product Type:</b>	Primary Antibodies
<b>Clone Name:</b>	YTS 156.7.7
<b>Applications:</b>	FC, IHC
<b>Reactivity:</b>	Mouse
<b>Host:</b>	Rat
<b>Isotype:</b>	IgG2b, kappa
<b>Clonality:</b>	Monoclonal
<b>Immunogen:</b>	Mouse CD8 (lyt-3).
<b>Specificity:</b>	The antibody recognises the B epitope of the mouse Lyt-3 (CD8) antigen. Has been used in CD8 subset depletion in vivo and flow cytometry studies.  This antibody recognizes murine CD8 beta.
<b>Formulation:</b>	PBS with 0.02% Proclin 300.
<b>Concentration:</b>	lot specific
<b>Conjugation:</b>	Unconjugated
<b>Storage:</b>	Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid freeze and thaw cycles.
<b>Stability:</b>	3 years from dispatch.
<b>Gene Name:</b>	CD8 antigen, beta chain 1
<b>Database Link:</b>	<a href="#">Entrez Gene 12526 Mouse P10300</a>
<b>Synonyms:</b>	Cd8b; Ly-3; Ly-C; Lyt-3



[View online »](#)