

Product datasheet for **TA386013**

VEGFA Rabbit Monoclonal Antibody [Clone ID: VG76e]

Product data:

Product Type:	Primary Antibodies
Clone Name:	VG76e
Applications:	ELISA, IHC, WB
Reactivity:	Human, Mouse, Primate
Host:	Rabbit
Isotype:	IgG, kappa
Clonality:	Monoclonal
Immunogen:	Human VEGF189 expressed in E. coli.
Specificity:	

VG76e will detect both native and denatured VEGF, an angiogenic growth factor and a prognostic indicator for cancer when detected in serum.

Formulation: PBS with 0.02% Proclin 300.

Concentration: lot specific

Conjugation: Unconjugated

Storage: Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid freeze and thaw cycles.

Stability: 3 years from dispatch.

Gene Name: vascular endothelial growth factor A

Database Link: [Entrez Gene 7422 Human P15692](#)

Synonyms: MGC70609; MVCD1; VEGF; VEGF-A; VPF

Note: This chimeric rabbit antibody was made using the variable domain sequences of the original Mouse IgG1 format, for improved compatibility with existing reagents, assays and techniques.



[View online »](#)