

Product datasheet for TA386031

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

TNFSF18 Rabbit Monoclonal Antibody [Clone ID: YGL386]

Product data:

Product Type: Primary Antibodies

Clone Name: YGL386

Applications: ELISA, FC

Reactivity: Mouse

Host: Rabbit

Isotype: IgG, kappa
Clonality: Monoclonal

Immunogen: Recombinant protein injected into DA rats.

Specificity: This antibody binds murine GITRL.

This antibody binds to GITRL and was used in the identification of GITRL and confirmation of

its function as the ligand for GITR.

Formulation: PBS with 0.02% Proclin 300.

Concentration: lot specific

Conjugation: Unconjugated

Storage: Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid

freeze and thaw cycles.

Stability: 3 years from dispatch.

Gene Name: tumor necrosis factor superfamily member 18

Database Link: Q9UNG2

Synonyms: AITRL; GITRL; hGITRL; MGC138237; TL6

Note: This chimeric rabbit antibody was made using the variable domain sequences of the original

Rat IgG1 format, for improved compatibility with existing reagents, assays and techniques.

