

Product datasheet for **TA386032**

TNFSF18 Rat Monoclonal Antibody [Clone ID: YGL386]

Product data:

| | |
|---------------|--|
| Product Type: | Primary Antibodies |
| Clone Name: | YGL386 |
| Applications: | ELISA, FC |
| Reactivity: | Mouse |
| Host: | Rat |
| Isotype: | IgG1, kappa |
| Clonality: | Monoclonal |
| Immunogen: | Recombinant protein injected into DA rats. |
| Specificity: | This antibody binds murine GITRL. |

This antibody binds to GITRL and was used in the identification of GITRL and confirmation of its function as the ligand for GTR.

| | |
|----------------|---|
| Formulation: | PBS with 0.02% Proclin 300. |
| Concentration: | lot specific |
| Conjugation: | Unconjugated |
| Storage: | Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid freeze and thaw cycles. |
| Stability: | 3 years from dispatch. |
| Gene Name: | tumor necrosis factor superfamily member 18 |
| Database Link: | Q9UNG2 |
| Synonyms: | AITRL; GITRL; hGITRL; MGC138237; TL6 |



[View online »](#)