

## Product datasheet for **TA386034**

### **TNFSF18 Rat Monoclonal Antibody [Clone ID: YGITR765]**

#### **Product data:**

<b>Product Type:</b>	Primary Antibodies
<b>Clone Name:</b>	YGITR765
<b>Applications:</b>	BI, FC
<b>Reactivity:</b>	Mouse
<b>Host:</b>	Rat
<b>Isotype:</b>	IgG2b, kappa
<b>Clonality:</b>	Monoclonal
<b>Immunogen:</b>	nB2/6TG cells transfected with mouse GITR injected into DA rats.
<b>Specificity:</b>	This antibody binds murine GITR.

A chimeric version of this antibody featuring a human IgG1 Fc-portion engineered to lack the N-linked glycosylation site at amino acid 297 was shown to enhance humoral immune responses by blocking GITR signallin (Ponte 2010).

<b>Formulation:</b>	PBS with 0.02% Proclin 300.
<b>Concentration:</b>	lot specific
<b>Conjugation:</b>	Unconjugated
<b>Storage:</b>	Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid freeze and thaw cycles.
<b>Stability:</b>	3 years from dispatch.
<b>Gene Name:</b>	tumor necrosis factor superfamily member 18
<b>Database Link:</b>	<a href="#">Q9UNG2</a>
<b>Synonyms:</b>	AITRL; GITRL; hGITRL; MGC138237; TL6



[View online »](#)