

## Product datasheet for **TA386035**

### Fibronectin (FN1) Human Monoclonal Antibody [Clone ID: A2]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	A2
Applications:	ELISA, IHC, WB
Reactivity:	Human
Host:	Human
Isotype:	IgG1, kappa
Clonality:	Monoclonal
Immunogen:	2D-PAGE gel isolated human Extra-Domain containing fibronectin from lysate of human melanoma COLO-38 cells.
Specificity:	This antibody bind to human Extra-Domain containing fibronectin.  Fibronectins are involved in cell adhesion, cell motility, opsonization, wound healing, and maintenance of cell shape. Extra domain-B containing fibronectin has been proposed marker of angiogenesis, has been shown to be expressed in a number of human cancers and in ocular neovascularization in patients with proliferative diabetic retinopathy.
Formulation:	PBS with 0.02% Proclin 300.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid freeze and thaw cycles.
Stability:	3 years from dispatch.
Gene Name:	fibronectin 1
Database Link:	<a href="#">Entrez Gene 2335 Human P02751</a>
Synonyms:	CIG; DKFZp686F10164; DKFZp686H0342; DKFZp686I1370; DKFZp686O13149; ED-B; FINC; FN; FNZ; GFND; GFND2; LETS; MSF; OTTHUMP00000206762; OTTHUMP00000206767



[View online »](#)

**Note:** This full-length, reformatted human antibody was made using the variable domain sequences of the original Human scFv format, for improved compatibility with existing reagents, assays and techniques.