

## Product datasheet for **TA386041**

### CD3E Mouse Monoclonal Antibody [Clone ID: 12F6]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	12F6
Applications:	FC
Reactivity:	Human
Host:	Mouse
Isotype:	IgG2a, kappa
Clonality:	Monoclonal
Immunogen:	Human CD3
Specificity:	This antibody binds to human CD3.

This antibody competes for CD3 binding with the therapeutic antibody OKT3 (Ab00122). 12F6 binding to CD3 on human T-cells stimulates the secretion of TNF-alpha, INF-gamma and IL-10. Humanized versions of this antibody are also available.

Formulation:	PBS with 0.02% Proclin 300.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid freeze and thaw cycles.
Stability:	3 years from dispatch.
Gene Name:	CD3e molecule
Database Link:	<a href="#">Entrez Gene 916 Human P07766</a>
Synonyms:	CD3-epsilon; FLJ18683; T3E; TCRE



[View online »](#)