

## Product datasheet for **TA386043**

### Integrin alpha 4 (ITGA4) Mouse Monoclonal Antibody [Clone ID: 21.6]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	21.6
Applications:	BI, ELISA
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1, kappa
Clonality:	Monoclonal
Immunogen:	Alpha-4 integrin.
Specificity:	This antibody binds to human alpha-4-integrin  This antibody binds to human alpha-4-integrin and can inhibit the binding of lymphocytes to endothelial cells, suppressing inflammation.
Formulation:	PBS with 0.02% Proclin 300.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid freeze and thaw cycles.
Stability:	3 years from dispatch.
Gene Name:	integrin subunit alpha 4
Database Link:	<a href="#">Entrez Gene 3676 Human P13612</a>
Synonyms:	CD49D; IA4; MGC90518; VLA-4
Note:	This reformatted mouse antibody was made using the variable domain sequences of the original Mouse IgG format, for improved compatibility with existing reagents, assays and techniques.



[View online »](#)