

Product datasheet for **TA386045**

BETVI Mouse Monoclonal Antibody [Clone ID: BV16]

Product data:

Product Type:	Primary Antibodies
Clone Name:	BV16
Applications:	ELISA, WB
Host:	Mouse
Isotype:	IgG1, kappa
Clonality:	Monoclonal
Immunogen:	Bet v 1-A, a major allergen of the European white birch.
Specificity:	This antibody binds to Bet v 1. This antibody binds to Bet v 1 via an epitope commonly used by allergenic IgE. Fab versions of this antibody have been used in co-crystallization studies of Bet v 1.
Formulation:	PBS with 0.02% Proclin 300.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid freeze and thaw cycles.
Stability:	3 years from dispatch.
Database Link:	P15494



[View online »](#)