

## Product datasheet for **TA386046**

### **BETVI Rabbit Monoclonal Antibody [Clone ID: BV16]**

#### **Product data:**

<b>Product Type:</b>	Primary Antibodies
<b>Clone Name:</b>	BV16
<b>Applications:</b>	ELISA, WB
<b>Host:</b>	Rabbit
<b>Isotype:</b>	IgG, kappa
<b>Clonality:</b>	Monoclonal
<b>Immunogen:</b>	Bet v 1-A, a major allergen of the European white birch.
<b>Specificity:</b>	This antibody binds to Bet v 1.  This antibody binds to Bet v 1 via an epitope commonly used by allergenic IgE. Fab versions of this antibody have been used in co-crystallization studies of Bet v 1.
<b>Formulation:</b>	PBS with 0.02% Proclin 300.
<b>Concentration:</b>	lot specific
<b>Conjugation:</b>	Unconjugated
<b>Storage:</b>	Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid freeze and thaw cycles.
<b>Stability:</b>	3 years from dispatch.
<b>Database Link:</b>	<a href="#">P15494</a>
<b>Note:</b>	This chimeric rabbit antibody was made using the variable domain sequences of the original Mouse IgG1 format, for improved compatibility with existing reagents, assays and techniques.



[View online »](#)