

## Product datasheet for **TA386048**

### **BETVI Rabbit Monoclonal Antibody [Clone ID: M0418]**

#### **Product data:**

<b>Product Type:</b>	Primary Antibodies
<b>Clone Name:</b>	M0418
<b>Applications:</b>	ELISA
<b>Host:</b>	Rabbit
<b>Clonality:</b>	Monoclonal
<b>Immunogen:</b>	Bet v 1-A, a major allergen of the European white birch.
<b>Specificity:</b>	<p>This antibody binds to Bet v 1, in particular to the naturally occurring isoforms Bet v 1.0101 and Bet v 1.0102, but not to Bet v 1.0112.</p> <p>This antibody binds to Bet v 1 via a characterized epitope and has been used in crystallographic studies.</p>
<b>Formulation:</b>	PBS with 0.02% Proclin 300.
<b>Concentration:</b>	lot specific
<b>Conjugation:</b>	Unconjugated
<b>Storage:</b>	Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid freeze and thaw cycles.
<b>Stability:</b>	3 years from dispatch.
<b>Database Link:</b>	<a href="#">P15494</a>
<b>Note:</b>	This full-length, chimeric rabbit antibody was made using the variable domain sequences of the original Human scFv format, for improved compatibility with existing reagents, assays and techniques.



[View online »](#)