

## Product datasheet for **TA386060**

### CD52 Mouse Monoclonal Antibody [Clone ID: CF1D12]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	CF1D12
Applications:	FC
Reactivity:	Human
Host:	Mouse
Isotype:	IgG3, kappa
Clonality:	Monoclonal
Specificity:	This antibody binds to human CD52.  This antibody binds to human CD52.
Formulation:	PBS with 0.02% Proclin 300.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid freeze and thaw cycles.
Stability:	3 years from dispatch.
Gene Name:	CD52 molecule
Database Link:	<a href="#">Entrez Gene 1043 Human P31358</a>
Synonyms:	CDw52; HE5



[View online »](#)