

Product datasheet for **TA386064**

Amyloid Precursor Protein (APP) Mouse Monoclonal Antibody [Clone ID: 60000000000]

Product data:

Product Type:	Primary Antibodies
Clone Name:	60000000000
Applications:	ELISA, EM, IHC, IP, WB
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1, kappa
Clonality:	Monoclonal
Immunogen:	Amyloid-beta peptide.
Specificity:	This antibody binds to amino acid residues 1-16 of beta amyloid. This antibody recognises human amyloid beta and can be used to detect amyloid in ELISAs, WB and tissue sections.
Formulation:	PBS with 0.02% Proclin 300.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid freeze and thaw cycles.
Stability:	3 years from dispatch.
Gene Name:	amyloid beta precursor protein
Database Link:	Entrez Gene 351 Human P05067
Synonyms:	A4; AAA; ABPP; AD1; APPI; CTFgamma; CVAP; OTTHUMP00000096096; PN-II; PN2; PreA4



[View online »](#)