

## Product datasheet for **TA386084**

### F Rabbit Monoclonal Antibody [Clone ID: RSHZ19 (Felvizumab)]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	RSHZ19 (Felvizumab)
Applications:	ELISA, FN
Reactivity:	Respiratory Syncytial Virus
Host:	Rabbit
Isotype:	IgG, kappa
Clonality:	Monoclonal
Immunogen:	This antibody was prepared by the humanization of an F protein-specific murine MAb, RSV19. More specifically, molecular techniques were used to insert the complementarity determining regions from murine RSV19 into a human IgG1 heavy and k light chain-variable domain framework.
Specificity:	This antibody is specific for the conserved fusion (F) protein of RSV.  Though shown to be safe and well-tolerated (Everitt, 1996) and protective against RSV infection in BALB/c mice and cotton rats (Wyde, 1995), this antibody did not demonstrate efficacy in a Phase 3 study (Reichert, 2008). In in vitro studies, palivizumab, another humanized Mab specific for the F protein of RSV, was 4 to 5-fold more potent than felvizumab in neutralising RSV (Scott, 1999).
Formulation:	PBS with 0.02% Proclin 300.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid freeze and thaw cycles.
Stability:	3 years from dispatch.
Note:	This chimeric rabbit antibody was made using the variable domain sequences of the original Human IgG1 format, for improved compatibility with existing reagents, assays and techniques.



[View online »](#)