

Product datasheet for TA386089

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

B7-1 (CD80) Human Monoclonal Antibody [Clone ID: IDEC-114 (Galiximab)]

Product data:

Product Type: Primary Antibodies

Clone Name: IDEC-114 (Galiximab)

Applications: BI, ELISA, FC, IHC, IP, WB

Reactivity: Human

Host: Human

Isotype: IgG1, lambda
Clonality: Monoclonal

Immunogen: Galiximab was prepared by immunizing cynomolgus monkeys with recombinant CD80

antigen. The variable regions of the light and heavy chains were then cloned by being incorporated into a cassette vector (N5LG1) containing human constant region genes and

subsequently transfected into the Dg44 CHO cell line (Hariharan K, 2013)

Specificity: Galiximab is a primatized mAb which consists of human constant and primate (cynomolgus

macaque) variable regions and binds specifically to CD80. CD80 is a surface glycoprotein and

a member of the B7 family of costimulatory molecules. CD80 antigen regulates T cell

activation (through interacting with CD28 or CD152) and is expressed transiently in antigenpresenting cells, T cells and normal B cells, and expressed consitutively on various subtypes

of B-cell lymphomas.

Galiximab acts to decrease in cell proliferation, inhibition of the constitutively active NF-κB pathway, increase in apoptosis and ADCC against various B-cell lymphoma cell lines. Galiximb affects CD28 and CD152 interactions with CD80. This antibody can also be used for FC, WB,

ELISA, IP and IHC.

Formulation: PBS with 0.02% Proclin 300.

Concentration: lot specific

Conjugation: Unconjugated

Storage: Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid

freeze and thaw cycles.

Stability: 3 years from dispatch.

Gene Name: CD80 molecule





B7-1 (CD80) Human Monoclonal Antibody [Clone ID: IDEC-114 (Galiximab)] - TA386089

Database Link: Entrez Gene 941 Human

P33681

Synonyms: B7; BB1; CD28LG; CD28LG1; LAB7