

## Product datasheet for **TA386096**

### Recombinant Anti-DDDDK-tag (Clone M2.1)

#### Product data:

Product Type:	Tag Antibodies
Product Name:	Recombinant Anti-DDDDK-tag (Clone M2.1)
Clone Name:	M2.1
Applications:	ELISA, IF, WB
Reactivity:	DDK (FLAG tag): DYKDDDDK
Host:	Rabbit
Clonality:	Monoclonal
Immunogen:	_x000D_ Interleukin 2 genetically fused to the peptided DYKDDDDK, corresponding to the DDDDK epitope tag.
Isotype:	IgG, kappa
Formulation:	PBS with 0.02% Proclin 300.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid freeze and thaw cycles.
Stability:	3 years from dispatch.
Note:	This chimeric rabbit antibody was made using the variable domain sequences of the original Mouse IgG1 format, for improved compatibility with existing reagents, assays and techniques.



[View online »](#)