

## Product datasheet for **TA386124**

### CD40L (CD40LG) Mouse Monoclonal Antibody [Clone ID: AT161-8]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	AT161-8
Applications:	FC, IF
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1, kappa
Clonality:	Monoclonal
Immunogen:	This antibody was raised by immunising mice with Chinese hamster ovary (CHO) cells expressing human CD40L.
Formulation:	PBS with 0.02% Proclin 300.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid freeze and thaw cycles.
Stability:	3 years from dispatch.
Gene Name:	CD40 ligand
Database Link:	<a href="#">Entrez Gene 959 Human P29965</a>
Synonyms:	CD40-L; CD40L; CD154; gp39; hCD40L; HIGM1; IGM; IMD3; T-BAM; TNFSF5; TRAP



[View online »](#)