

Product datasheet for **TA386129**

Recombinant Anti-Podoplanin (MAP tag) (Clone PMab-1)

Product data:

Product Type:	Tag Antibodies
Product Name:	Recombinant Anti-Podoplanin (MAP tag) (Clone PMab-1)
Clone Name:	PMab-1
Applications:	BI, ELISA, IHC, WB
Reactivity:	MAP tag
Host:	Rat
Clonality:	Monoclonal
Immunogen:	This antibody was raised by immunising rats with 14-residue synthetic peptide mpp3851, which corresponds to amino acids 38-51 (GDGMVPPGIEDKIT) of the platelet-aggregation-stimulating (PLAG) domain of mouse podoplanin. Spleen cells were then harvested, and fused to P3U1 cells to generate stable hybridomas.
Isotype:	IgG2a, kappa
Formulation:	PBS with 0.02% Proclin 300.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid freeze and thaw cycles.
Stability:	3 years from dispatch.
Gene Name:	podoplanin



[View online »](#)