

Product datasheet for **TA386139**

Tnfsf9 Rat Monoclonal Antibody [Clone ID: TKS-1]

Product data:

Product Type:	Primary Antibodies
Clone Name:	TKS-1
Applications:	BI, FC
Reactivity:	Mouse
Host:	Rat
Isotype:	IgG2a, kappa
Clonality:	Monoclonal
Immunogen:	Murine B-cell lymphoma 2PK-3 cells injected into SD rats.
Specificity:	This antibody binds to murine 4-1BBL and competes with binding of 4-1BB. This antibody has been used to detect 4-1BBL by flow cytometry. It is able to block 4-1BBL interacting with 4-1BB and prevent signalling.
Formulation:	PBS with 0.02% Proclin 300.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid freeze and thaw cycles.
Stability:	3 years from dispatch.
Gene Name:	tumor necrosis factor (ligand) superfamily, member 9
Database Link:	Entrez Gene 21950 Mouse P41274
Synonyms:	4-1BB-L; 4-1BBL; CD137L



[View online »](#)