

Product datasheet for **TA386141**

Icosl Rat Monoclonal Antibody [Clone ID: HK5.3]

Product data:

Product Type:	Primary Antibodies
Clone Name:	HK5.3
Applications:	BI, FC
Reactivity:	Mouse
Host:	Rat
Isotype:	IgG2a, kappa
Clonality:	Monoclonal
Immunogen:	This antibody was raised by immunizing SD rat with mouse B7h-transfected L cells and fusing immune splenocytes with P3U1 myeloma cells and screened for binding to mouse B7h-transfected NRK cells.
Specificity:	<p>This antibody specifically binds to mouse ICOSL and competes with binding of ICOS.</p> <p>This antibody has been used to detect murine ICOSL by flow cytometry. It is able to block ICOSL interacting with ICOSL and inhibit T-cell proliferation in a collagen type II (CII)-induced arthritis model.</p>
Formulation:	PBS with 0.02% Proclin 300.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid freeze and thaw cycles.
Stability:	3 years from dispatch.
Gene Name:	icos ligand
Database Link:	Entrez Gene 50723 Mouse Q9JHJ8
Synonyms:	AU044799; B7-H2; B7h; B7RP-1; BG071784; GI50; GL50; GL50-B; ICOS-L; Icoslg; LICOS; Ly115I



[View online »](#)