

Product datasheet for **TA386202**

Cd80 Rabbit Monoclonal Antibody [Clone ID: RM80]

Product data:

Product Type:	Primary Antibodies
Clone Name:	RM80
Applications:	BI, FC, FN, IHC, WB
Reactivity:	Mouse
Host:	Rabbit
Isotype:	IgG, kappa
Clonality:	Monoclonal
Immunogen:	This antibody was raised by immunising Sprague Dawley rats with a mouse B cell line, BCL1 and fusing immune splenocytes with P3U1 myeloma cells.



[View online »](#)

Specificity:	<p>This antibody is specific for murine CD80, a ~55 kDa member of the Ig superfamily. CD80 is expressed by macrophages, dendritic cells, activated B cells, and activated T cells. CD80 has high affinity for binding to two T cell surface antigens, CD28 and CD152 (CTLA-4). The interaction of CD28 and CD152 with CD80 is crucial in T-B cell communication leading to the activation of T and B cells, respectively.</p> <p>This antibody has been used in multiple FACS experiments for various immunological applications, such as to highlight the robustness of the T cell priming process by activated Nlr5-deficient dendritic cells (Rota et al, 2016), to demonstrate that allogeneic IgG combined with dendritic cell stimuli could induce anti-tumor T cell immunity (Carmi et al, 2015), to confirm the role of aquaporin-9-expressing neutrophils for the development of contact hypersensitivity (Moniaga et al, 2015), and to support the therapeutic potential of targeting B-cell Syk signaling in chronic graft-versus-host disease (Flynn et al, 2015). This antibody has also been used in Western Blot to study the effect of interleukin 4-induced gene-1 (IL411) overexpression on the expression of M2 markers and M1-associated cytokines (Yue et al, 2015), as well as in in vitro blocking experiments to confirm the role of CD28 signals for IL-2-mediated survival and proliferation of both CD27(+) and CD27(-) $\gamma\delta$ T cell subsets (Ribot et al, 2012). Furthermore, this antibody has been used, together with the anti-CD86 antibody, in numerous in vivo studies, such as in a NZB/W F1 mouse lupus model to clarify the role of CD80 and CD86 in murine Ab-mediated autoimmunity (Nakajima et al, 1995), in graft bearing recipient BALB/c mice to evaluate the roles of CD80 and CD86 costimulatory molecules in corneal allograft rejection (Kagaya et al, 2002), and in transgenic BALB/c mice to demonstrate the critical contribution of CD80 and CD86 in anterior chamber-associated immune deviation (Tsukahara et al, 2005).</p>
Formulation:	PBS with 0.02% Proclin 300.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid freeze and thaw cycles.
Stability:	3 years from dispatch.
Gene Name:	CD80 antigen
Database Link:	Entrez Gene 12519 Mouse Q00609
Synonyms:	B7; BB1; CD28LG; CD28LG1; LAB7
Note:	This chimeric rabbit antibody was made using the variable domain sequences of the original Rat IgG2a format, for improved compatibility with existing reagents, assays and techniques.