

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA386204

Cd86 Rabbit Monoclonal Antibody [Clone ID: PO.3]

Product data:

Product Type:	Primary Antibodies
Clone Name:	PO.3
Applications:	BI, FC, FN, IHC, WB
Reactivity:	Mouse
Host:	Rabbit
lsotype:	IgG, kappa
Clonality:	Monoclonal
Immunogen:	This antibody was raised by immunising Sprague Dawley rats with a mouse B cell line, BCL1 and fusing immune splenocytes with P3U1 myeloma cells.
Specificity:	This antibody is specific for murine CD86, an ~80 kDa surface receptor of the B7 family of costimulatory molecules. In addition to CD80, CD86 is a counter-receptor for the T cell surface molecules CD28 and CD152 (CTLA-4). This interaction plays a critical role in T-B crosstalk, T cell costimulation, autoantibody production and Th2-mediated Ig production. This antibody has been used in multiple FACS experiments for various immunological applications, such as to demonstrate how the C-type lectin receptor DCIR modulates immunity to tuberculosis (Troegeler et al, 2017), to elucidate the mechanism of Th17 cell differentiation by cholera toxin (Kang et al, 2016), and to evaluate the role of the spleen tyrosine kinase in graft-versus-host disease (Flynn et al, 2015). This antibody has also been used in immunofluorescence to investigate the potential of peritoneal cavity B cells to alleviate murine colitis (Maseda et al, 2013), and in Western Blot to study the effect of interleukin 4-induced gene-1 (IL4I1) overexpression on the expression of M2 markers and M1-associated cytokines (Yue et al, 2015). Furthermore, this antibody has been used, together with the anti-CD80 antibody, in numerous in vivo studies, such as in a NZB/W F1 mouse lupus model to clarify the role of CD80 and CD86 in murine Ab-mediated autoimmunity (Nakajima et al, 1995), in graft bearing recipient BALB/c mice to evaluate the roles of CD80 and CD86 in anterior
	chamber-associated immune deviation (Tsukahara et al, 2005).
Formulation:	PBS with 0.02% Proclin 300.



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

Cd86 Rabbit Monoclonal Antibody [Clone ID: PO.3] – TA386204

Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid freeze and thaw cycles.
Stability:	3 years from dispatch.
Gene Name:	CD86 antigen
Database Link:	<u>Entrez Gene 12524 Mouse</u> <u>P42082</u>
Synonyms:	B7-2; B70; BU63; CD28LG2; FUN-1; LAB72; MGC34413
Note:	This chimeric rabbit antibody was made using the variable domain sequences of the original Rat IgG2b format, for improved compatibility with existing reagents, assays and techniques.

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US