

Product datasheet for TA386228

Timd4 Rabbit Monoclonal Antibody [Clone ID: 5G3]

5G3

Mouse

Rabbit

lgG, kappa

Monoclonal

Primary Antibodies

BI, ELISA, FC, IF, IHC, WB

4:lg fusion protein in complete Freund's adjuvant.

This antibody was raised by immunising female Lewis rats (Harlan Sprague-Dawley) with Tim-

Product data:

Product Type:

Clone Name:

Applications:

Reactivity:

Host:

Isotype:

Clonality:

Immunogen:

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

| | Timd4 Rabbit Monoclonal Antibody [Clone ID: 5G3] – TA386228 |
|----------------|---|
| Specificity: | This antibody is specific for murine TIM-4, which unlike other TIM family members, is not expressed on T lymphocytes but rather found on antigen presenting cells such as lymphoid and myeloid DC's and macrophages. TIM-4 has been shown to interact with TIM-1 thereby suggesting a role of TIM-4 in the balance of Th1/Th2 regulation. Additionally TIM-4 has been shown to be a receptor for phosphatidylserine and may play a role in engulfment of apoptotic cells. |
| | This antibody (Clone 5G3), together with two others (clones 3A1 and 3H11), was originally produced from the subcloned hybridomas; all of these 3 produced mAbs bound specifically to murine TIM-4 and not to any other TIM molecule, as demonstrated by both FACS and ELISA analyses (Rodriguez-Manzanet et al, 2008). This antibody has been used in various FACS analyses, for example, to show that TIM-1 and TIM-4 glycoproteins bind phosphatidylserine and mediate uptake of apoptotic cells (Kobayashi et al, 2007), to underscore the requirement of IL-10 from marginal zone precursor B cell subset for costimulatory blockade-induced transplantation tolerance (Lal et al, 2015), and to uncover a mechanism by which the body adapts to fluctuations in erythrocyte integrity (Theurl et al, 2016). This antibody has also been used in immunohistochemistry to discover a new tissue macrophage subset in the thymus that is dependent on transcription factor NR4A1 for development (Tacke et al, 2015), as well as in immunofluorescence to elucidate the mechanisms that control the heterogeneity of tissue macrophages (A-Gonzalez et al, 2013). In addition, the anti-tumour effects of anti-TIM-4 antibodies have been demonstrated in various murine melanoma and colon adenocarcinoma models (Baghdadi et al, Cancer Immunol. Immunother. 2013; Baghdadi et al, Immunity. 2013). |
| Formulation: | PBS with 0.02% Proclin 300. |
| Concentration: | lot specific |
| Conjugation: | Unconjugated |
| Storage: | Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid freeze and thaw cycles. |
| Stability: | 3 years from dispatch. |
| Gene Name: | T cell immunoglobulin and mucin domain containing 4 |
| Database Link: | <u>Entrez Gene 276891 Mouse</u> <u>Q6U7R4</u> |
| Synonyms: | FLJ27515; SMUCKLER; TIM-4; TIM4; TIMD-4 |
| Note: | This chimeric rabbit antibody was made using the variable domain sequences of the original Rat lgG1 format, for improved compatibility with existing reagents, assays and techniques. |

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US