

## Product datasheet for **TA386239**

### **Tnfsf8 Rat Monoclonal Antibody [Clone ID: RM153]**

#### **Product data:**

**Product Type:** Primary Antibodies

**Clone Name:** RM153

**Applications:** BI, FC, FN

**Reactivity:** Mouse

**Host:** Rat

**Isotype:** IgG2b, kappa

**Clonality:** Monoclonal

**Immunogen:** This antibody was raised by immunising Sprague Dawley rats with the CHO cells stably transfected with murine CD153 and fusing the splenocytes with P3U1 myeloma cells.

**Specificity:** This antibody is specific for murine CD153, a 40-kDa type II transmembrane protein member of the TNF superfamily. Interaction of CD30 with CD153 plays an important role in T cell activation/costimulation.

This antibody has been used in various FACS analyses, for instance, to demonstrate that lymphoid tissue inducer-like cells are an innate source of IL-17 and IL-22 (Takatori et al, 2009), and to suggest that CD30 ligand/CD30 plays a critical role in Th17 differentiation in mice (Sun et al, 2010). Importantly, the pharmacological effects of this rat anti-mouse CD153 mAb (clone RM153; IgG2b isotype) have been reported in a number of settings. For example, this antibody has been shown to prevent or delay the spontaneous development of diabetes in non-obese diabetic mice when treatment is initiated in young mice (Chakrabarty et al, 2003), and to prolong the survival of mice in a model of CD4+ T-cell-mediated lethal graft-versus-host disease induced by injection of C57BL/6 cells into MHC class II-only disparate bm12 recipients (Blazar et al, 2004). This antibody has also been used as parts of in vivo functional assays to demonstrate that early CD30 signaling is critical for Treg-mediated acute graft-versus-host disease protection after major MHC-mismatch bone marrow transplantation (Zeiser et al, 2007), and to provide direct evidence that pathogenic memory T cells are amenable to suppression in an antigen-specific manner while identifying CD30 as a molecule that is critical for the regulation of memory T cell responses (Dai et al, 2004).

**Formulation:** PBS with 0.02% Proclin 300.

**Concentration:** lot specific



[View online »](#)

<b>Conjugation:</b>	Unconjugated
<b>Storage:</b>	Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid freeze and thaw cycles.
<b>Stability:</b>	3 years from dispatch.
<b>Gene Name:</b>	tumor necrosis factor (ligand) superfamily, member 8
<b>Database Link:</b>	<a href="#">Entrez Gene 21949 Mouse P32972</a>
<b>Synonyms:</b>	CD30-L; CD30L; CD30LG; CD153; MGC138144