

## Product datasheet for **TA386245**

### Ifngr1 Rat Monoclonal Antibody [Clone ID: GR20]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	GR20
Applications:	BI, FC, Neutralize
Reactivity:	Mouse
Host:	Rat
Isotype:	IgG2a, kappa
Clonality:	Monoclonal
Immunogen:	This antibody was raised by immunising rats with interferon gamma receptor.
Specificity:	<p>This antibody recognises IFN-<math>\gamma</math>R<math>\alpha</math>, which is a class II cytokine receptor family member that serves as a IFN-<math>\gamma</math>-binding chain associated with the IFN-<math>\gamma</math> <math>\beta</math>.</p> <p>This antibody is recommended as an efficient blocking agent for the interferon gamma receptor (Basu et al. 1988). It has been also used in the flow cytometric analyses of this receptor protein.</p>
Formulation:	PBS with 0.02% Proclin 300.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid freeze and thaw cycles.
Stability:	3 years from dispatch.
Gene Name:	interferon gamma receptor 1
Database Link:	<a href="#">Entrez Gene 15979 Mouse P15261</a>
Synonyms:	CD119; CDw119; FLJ45734; IFN-gamma-R1; IFNGR



[View online »](#)