

## **Product datasheet for TA386258**

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

# **Tnfsf11 Rabbit Monoclonal Antibody [Clone ID: IK22-5]**

### **Product data:**

**Product Type:** Primary Antibodies

Clone Name: IK22-5

**Applications:** FC, IHC, Neutralize

**Reactivity:** Mouse **Host:** Rabbit

Isotype:IgG, kappaClonality:Monoclonal

Immunogen: Thus antibody was raised by immunising an SD rat with murine RANKL-transfected baby

hamster kidney cells.

**Specificity:** This antibody is specific for murine RANKL.

This antibody has been used in immunohistochemical analysis of frozen Peyer's patch tissue sections (Kobayashi et al, 2013). It has also been used in flow cytometric analysis of RANKL cell surface expression by Th1 and Th17 cells (Nakae et al, 2007), RANKL-transfected human embryonic kidney cells (Miyahira et al, 2003) and murine T hybridoma 2B4.11 cells (Fionda et al, 2007). Additionally, this RANKL-neutralising antibody has been used to assess the effects of RANKL inhibition on myoepithelial cell proliferation in mammary organoid culture (Haslam et al, 2008). Treatment of collagen-induced arthritis (CIA) mice with this antibody has been

shown to reduce bone loss at the site of inflammation (Kamijo et al, 2006), while

administration to C57BL/6J mice induced transient tibia and skull bone growth arrest and

irreversible blockage of tooth eruption (Lezot et al, 2015).

**Formulation:** PBS with 0.02% Proclin 300.

**Concentration:** lot specific

**Conjugation:** Unconjugated

**Storage:** Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid

freeze and thaw cycles.

**Stability:** 3 years from dispatch.

**Gene Name:** tumor necrosis factor (ligand) superfamily, member 11





### Tnfsf11 Rabbit Monoclonal Antibody [Clone ID: IK22-5] - TA386258

Database Link: Entrez Gene 21943 Mouse

<u>O35235</u>

Synonyms: CD254; hRANKL2; ODF; OPGL; OPTB2; OTTHUMP00000178585; RANKL; sOdf; TRANCE

**Note:** This chimeric rabbit antibody was made using the variable domain sequences of the original

rat IgG2a format, for improved compatibility with existing reagents, assays and techniques.