

Product datasheet for **TA386287**

IL9 Rat Monoclonal Antibody [Clone ID: RM9A4]

Product data:

Product Type:	Primary Antibodies
Clone Name:	RM9A4
Applications:	ELISA, FC
Reactivity:	Mouse
Host:	Rat
Isotype:	IgG1, kappa
Clonality:	Monoclonal
Immunogen:	This antibody was raised by immunising rats with recombinant full-length murine IL-9.
Specificity:	This antibody is specific for murine IL-9.

This antibody, when fluorochrome-conjugated, has been used to achieve intracellular IL-9 staining of T cells, analysed using flow cytometry (Chang et al, 2010; Kim et al, 2010; Horka et al, 2012; Gomez-Rodriguez et al, 2016). This antibody has also been used in a sandwich ELISA to measure IL-9 secretion in supernatants collected from restimulated T cells (Jager et al, 2009).

Formulation:	PBS with 0.02% Proclin 300.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid freeze and thaw cycles.
Stability:	3 years from dispatch.
Gene Name:	interleukin 9
Database Link:	Entrez Gene 16198 Mouse P15247
Synonyms:	HP40; IL-9; P40



[View online »](#)