

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA386296

Lag3 Rabbit Monoclonal Antibody [Clone ID: C9B7W]

Product data:

Product Type:	Primary Antibodies
Clone Name:	C9B7W
Applications:	BI, FC, IP
Reactivity:	Mouse
Host:	Rabbit
lsotype:	lgG, kappa
Clonality:	Monoclonal
Immunogen:	This antibody was raised by immunising Lewis rats with a murine LAG-3:Ig fusion protein.
Specificity:	This antibody is specific for an epitope in the D2 domain of murine LAG-3.
	This antibody has been used in flow cytometric analysis of splenocytes (Workman et al, 2005), T cell hybridomas expressing wild-type or chimeric CD4:murine LAG-3 (Workman et al, 2002) and CD4+ regulatory T cells (Chien et al, 2017). It has also been used in immunoprecipitation analysis of murine T hybridoma 3A9 cells transduced with murine LAG-3 (Li et al, 2004). This antibody blocks LAG-3 function in vitro, inhibiting antigen-induced IL-2 production in LAG-3- expressing T cell hybridomas, but does not block LAG-3:MHC class II interaction (Workman et al, 2002). Treatment with this antibody also blocks LAG-3 function in vivo, leading to accelerated homeostatic expansion of WT T cells (Workman et al, 2005). Administration of this LAG-3-blocking antibody, in combination with PD-L1 blockade, synergistically improves the survival of lymphodepleted myeloma-bearing mice (Jing et al, 2015).
Formulation:	PBS with 0.02% Proclin 300.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid freeze and thaw cycles.
Stability:	3 years from dispatch.
Gene Name:	lymphocyte-activation gene 3
Database Link:	<u>Entrez Gene 16768 Mouse</u> <u>Q61790</u>



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

	Lag3 Rabbit Monoclonal Antibody [Clone ID: C9B7W] – TA386296
Synonyms:	CD223; FDC; LAG-3
Note:	This chimeric rabbit antibody was made using the variable domain sequences of the original Rat IgG1 format, for improved compatibility with existing reagents, assays and techniques.

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US