

## Product datasheet for **TA386305**

### IL2rg Rat Monoclonal Antibody [Clone ID: 3000000000000]

#### Product data:

<b>Product Type:</b>	Primary Antibodies
<b>Clone Name:</b>	3000000000000
<b>Applications:</b>	BI, FC, IP
<b>Reactivity:</b>	Mouse
<b>Host:</b>	Rat
<b>Isotype:</b>	IgG2b, kappa
<b>Clonality:</b>	Monoclonal
<b>Immunogen:</b>	This antibody was raised by immunising Lewis rats with transfected rat myeloma cell line YB2/0, expressing murine cytoplasmic-tailless IL-2R common gamma-chain (YB2/0-myt).
<b>Specificity:</b>	<p>This conformation-dependent antibody binds to the murine common gamma-chain, shared by the IL-2, IL-4, IL-7, IL-9, IL-15, and IL-21 receptor complexes. It recognises an epitope located in a membrane-proximal portion of the common gamma-chain, near Arg-197 of domain 2. This antibody does not cross-react with COS7 cells transfected with rat common gamma-chain cDNA (Olosz et al, 2002).</p> <p>This antibody has been used in flow cytometric analysis of common gamma-chain surface expression by IL-7-stimulated T cells (Goh et al, 2017) and COS7 monkey kidney cells (Olosz et al, 2002). This antibody has also been used to immunoprecipitate murine common gamma-chain-containing complexes from COS7 cell extracts (Olosz et al, 2002). Additionally, as this antibody is conformationally sensitive, it has been used to confirm the correct folding of soluble murine common gamma chain mutants (Meißner et al, 2001). This antagonistic antibody is most effective at inhibiting IL-4 bioactivity (Malek et al, 1998). Treatment of recipient mice with this antibody, along with another common gamma-chain-blocking MAb (4G3), significantly prolongs islet allograft survival (Li et al, 2000). Indeed, blocking the common gamma chain with this antibody has been shown to induce apoptotic cell death of BAF-B03 cells and primary anti-CD3 activated T lymphoblasts in vitro, and of alloreactive T cells in vivo (Li et al, 2000).</p>
<b>Formulation:</b>	PBS with 0.02% Proclin 300.
<b>Concentration:</b>	lot specific
<b>Conjugation:</b>	Unconjugated



[View online »](#)

<b>Storage:</b>	Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid freeze and thaw cycles.
<b>Stability:</b>	3 years from dispatch.
<b>Gene Name:</b>	interleukin 2 receptor, gamma chain
<b>Database Link:</b>	<a href="#">Entrez Gene 16186 Mouse P34902</a>
<b>Synonyms:</b>	CD132; Gamma-C; IMD4; p64; SCIDX; SCIDX1