

## Product datasheet for **TA386320**

### FCER1G Mouse Monoclonal Antibody [Clone ID: IE7]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	IE7
Applications:	IF
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1, kappa
Clonality:	Monoclonal
Immunogen:	Recombinant human protein, N-terminus.
Specificity:	This antibody binds to human Fc epsilon receptor, but does not compete with IgE binding.  IE7 binds to human FcεR1a, but does not compete with IgE binding or binding of the antibody IB10 (Ab01396).
Formulation:	PBS with 0.02% Proclin 300.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid freeze and thaw cycles.
Stability:	3 years from dispatch.
Gene Name:	Fc fragment of IgE receptor Ig
Database Link:	<a href="#">Entrez Gene 2207 Human P30273</a>
Synonyms:	FCRG



[View online »](#)