

Product datasheet for **TA386356**

CD41 (ITGA2B) Mouse Monoclonal Antibody [Clone ID: BAG 75E5]

Product data:

Product Type:	Primary Antibodies
Clone Name:	BAG 75E5
Applications:	IHC
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1, kappa
Clonality:	Monoclonal
Immunogen:	This antibody was raised by immunising mice with human CD41.
Specificity:	This antibody recognises CD41 antigen. This antibody is recommended for the detection and identification of megakaryocytes and thrombocytes. It was used to analyse megakaryocytes in the fetal hepatoblastoma tissue samples (von Schweinitz et al., 1995).
Formulation:	PBS with 0.02% Proclin 300.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid freeze and thaw cycles.
Stability:	3 years from dispatch.
Gene Name:	integrin subunit alpha 2b
Database Link:	Entrez Gene 3674 Human P08514
Synonyms:	CD41; CD41B; GP2B; GPIIb; GTA; HPA3; ITGAB



[View online »](#)