

Product datasheet for **TA386464**

PRDM4 Human Monoclonal Antibody [Clone ID: RAB-C367]

Product data:

Product Type:	Primary Antibodies
Clone Name:	RAB-C367
Applications:	ChIP, ELISA, IF
Reactivity:	Human
Host:	Human
Isotype:	Fab fragment, kappa
Clonality:	Monoclonal
Immunogen:	This antibody was obtained by recombinant antibody (rAb) phage display recognizing PRDM4 protein under non-denaturing conditions.
Specificity:	<p>This antibody recognizes PRDM4 (PR domain zinc finger protein 4). It binds to a folded domain, amino acids 390-540. PRDM4 is a transcription factor involved in cell differentiation.</p> <p>This antibody has been used and recommended for the following methods: ELISA; IF; IP-MS; ChIP.</p>
Formulation:	PBS with 0.02% Proclin 300.
Concentration:	lot specific
Conjugation:	His Tag
Storage:	Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid freeze and thaw cycles.
Stability:	3 years from dispatch.
Gene Name:	PR domain 4
Database Link:	Entrez Gene 11108 Human Q9UKN5
Synonyms:	MGC45046; PFM1



[View online »](#)