

Product datasheet for **TA386468**

PRMT3 Human Monoclonal Antibody [Clone ID: RAB-C372]

Product data:

Product Type:	Primary Antibodies
Clone Name:	RAB-C372
Applications:	ChIP, ELISA
Reactivity:	Human
Host:	Human
Isotype:	Fab fragment, kappa
Clonality:	Monoclonal
Immunogen:	This antibody was obtained by recombinant antibody (rAb) phage display recognizing PRMT3 protein under non-denaturing conditions.
Specificity:	<p>This antibody recognizes PRMT3 (Protein arginine N-methyltransferase 3). It binds to a folded domain, amino acids 211-531. PRMT3 methylates (mono and asymmetric dimethylation) the guanidino nitrogens of arginyl residues in some proteins.</p> <p>This antibody has been used and recommended for the following methods: ELISA; IP-MS; ChIP.</p>
Formulation:	PBS with 0.02% Proclin 300.
Concentration:	lot specific
Conjugation:	His Tag
Storage:	Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid freeze and thaw cycles.
Stability:	3 years from dispatch.
Gene Name:	protein arginine methyltransferase 3
Database Link:	Entrez Gene 10196 Human O60678
Synonyms:	HRMT1L3



[View online »](#)