

## Product datasheet for **TA386478**

### FOXP4 Human Monoclonal Antibody [Clone ID: RAB-S40]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	RAB-S40
Applications:	ELISA, FC, IF
Reactivity:	Human
Host:	Human
Isotype:	Fab fragment, kappa
Clonality:	Monoclonal
Immunogen:	This antibody was obtained by recombinant antibody (rAb) phage display recognizing FOXP4 protein under non-denaturing conditions.
Specificity:	This antibody recognizes FOXP4 (Forkhead box protein P4). It binds to a folded domain, amino acids 302-541. FOXP4 is a transcriptional repressor that represses lung-specific expression.  This antibody has been used and recommended for the following methods: ELISA; spiked-IP; IF; FC.
Formulation:	PBS with 0.02% Proclin 300.
Concentration:	lot specific
Conjugation:	His Tag
Storage:	Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid freeze and thaw cycles.
Stability:	3 years from dispatch.
Gene Name:	forkhead box P4
Database Link:	<a href="#">Entrez Gene 116113 Human Q8IVH2</a>
Synonyms:	FKHLA; FLJ40908; FLJ44184; hFKHLA; OTTHUMP00000043536



[View online »](#)