

Product datasheet for **TA386485**

CD38 Rabbit Monoclonal Antibody [Clone ID: OKT10]

Product data:

Product Type:	Primary Antibodies
Clone Name:	OKT10
Applications:	FC, IP
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG, kappa
Clonality:	Monoclonal
Immunogen:	The original antibody was generated by immunizing CAF1 mice with human thymocytes.
Specificity:	This antibody recognizes and antigenic structure of 45,000 daltons type II transmembrane glycoprotein present on human precursor B cells and T cells, activated T cells, natural killer cells and some monocytes, myeloblasts and erythroblasts. It is also reported to strongly react with surface membrane of plasma cells.

This antibody reacts with various lymphoid cell populations and thymocytes. OKT10 along with other antibodies was used to describe the stages of intrathymic development of human thymocytes. It was reported that all the thymocytes initially bear the OKT10 marker. As the thymocytes are exported to peripheral T-cell compartment they lose the OKT10 marker, as OKT10 could not react with peripheral T-cells (PMID:6966400). OKT10 is reported to bind to leukaemias and normal marrow cells of the myeloid and B cell lineages (PMID:6169674). OKT10 was used to study final B cell maturation since the OKT10 antigen is expressed on mature lymphocytes and lymphoplasmacytic cells. OKT10 reactivity is inversely proportional to the presence of Ia-like antigen expression. Therefore, these two antigens can be employed as complementary opposite reactants for the study of mature B cell malignancies, such as CLL, multiple myeloma, and Waldenström malignancy. (PMID:6813349). This antibody was used for determination of highly purified human NK cells by flow cytometry (PMID:6975321). This antibody was used for the identification of a CD38-like 78-Kilodalton soluble protein released from B Cell Lines derived from patients with X-linked Agammaglobulinemia (PMID:9637716)

Formulation: PBS with 0.02% Proclin 300.

Concentration: lot specific



[View online »](#)

Conjugation:	Unconjugated
Storage:	Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid freeze and thaw cycles.
Stability:	3 years from dispatch.
Gene Name:	CD38 molecule
Database Link:	Entrez Gene 952 Human P28907
Synonyms:	T10
Note:	This chimeric rabbit antibody was made using the variable domain sequences of the original Mouse IgG1 format, for improved compatibility with existing reagents, assays and techniques.