

Product datasheet for **TA386626M**

ZMYND10 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200
Reactivity:	Mouse, Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant Human Zinc finger MYND domain-containing protein 10 protein (201-440AA)
Formulation:	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Concentration:	lot specific
Purification:	Antigen Affinity Purified
Conjugation:	Unconjugated
Storage:	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Stability:	1 year from dispatch.
Database Link:	O75800
Background:	Required for motile ciliary function. Probably involved in axonemal assembly of inner and outer dynein arms (IDA and ODA, respectively) for proper axoneme building for cilia motility. May act by indirectly regulating transcription of dynein proteins.




[View online »](#)


Product images:



Immunohistochemistry of paraffin-embedded human testis tissue using [TA386626] at dilution of 1:100



Immunohistochemistry of paraffin-embedded human lung cancer using [TA386626] at dilution of 1:100



Western blot

All lanes: ZMYND10 antibody at 8 μ g/ml

Lane 1: A549 whole cell lysate

Lane 2: Mouse kidney tissue

Lane 3: Mouse skeletal muscle tissue

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 51, 50 kDa

Observed band size: 51 kDa