

Product datasheet for **TA386627**


ANXA9 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	Recommended dilution: WB:1:500–1:2000, IHC:1:20–1:200
Reactivity:	Mouse, Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant Human Annexin A9 protein (8–345AA)
Formulation:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Concentration:	lot specific
Purification:	>95%, Protein G purified
Conjugation:	Unconjugated
Storage:	Upon receipt, store at –20°C or –80°C. Avoid repeated freeze.
Stability:	1 year from dispatch.
Database Link:	<u>O76027</u>
Background:	Low affinity receptor for acetylcholine known to be targeted by disease-causing pemphigus vulgaris antibodies in keratinocytes.



Product images:



Western blot

All lanes: ANXA9 antibody at 8 μ g/ml

Lane 1: Mouse kidney tissue


Lane 2: A431 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution


Predicted band size: 39 kDa

Observed band size: 39 kDa



Immunohistochemistry of paraffin-

embedded human placenta tissue using
TA386627 at dilution of 1:100



Immunohistochemistry of paraffin-embedded human thyroid tissue using TA386627 at dilution of 1:100