

Product datasheet for **TA386691M**


GPX4 Rabbit Polyclonal Antibody

Product data:


Product Type:	Primary Antibodies
Applications:	IF, IHC, WB
Recommended Dilution:	Recommended dilution: WB:1:500-1:5000, IHC:1:20-1:200, IF:1:50-1:200
Reactivity:	Mouse, Rat, Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant Human Phospholipid hydroperoxide glutathione peroxidase, mitochondrial protein (74-197AA)
Formulation:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Concentration:	lot specific
Purification:	>95%, Protein G purified
Conjugation:	Unconjugated
Storage:	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Stability:	1 year from dispatch.
Database Link:	P36969
Background:	Protects cells against membrane lipid peroxidation and cell death. Required for normal sperm development and male fertility. Could play a major role in protecting mammals from the toxicity of ingested lipid hydroperoxides. Essential for embryonic development. Protects from radiation and oxidative damage. Essential for maturation and survival of photoreceptor cells. Plays a role in a primary T cell response to viral and parasitic infection by protecting T cells from ferroptosis, a cell death resulting from an iron-dependent accumulation of lipid reactive oxygen species, and by supporting T cell expansion.

[View online »](#)


Product images:



Immunohistochemistry of paraffin-embedded human kidney tissue using [TA386691] at dilution of 1:100



Immunofluorescent analysis of A549 cells using [TA386691] at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)



Western Blot

Positive WB detected in: A549 whole cell lysate, MCF-7 whole cell lysate, Mouse kidney tissue, Rat heart tissue

All lanes: GPX4 antibody at 3.5µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 23, 20 kDa

Observed band size: 23 kDa