

Product datasheet for **TA386732M**

TOB2 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Reactivity:	Mouse, Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant Human Protein Tob2 protein (1-70AA)
Formulation:	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Concentration:	lot specific
Purification:	Antigen Affinity Purified
Conjugation:	Unconjugated
Storage:	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Stability:	1 year from dispatch.
Database Link:	Q14106
Background:	Anti-proliferative protein inhibits cell cycle progression from the G0/G1 to S phases.




[View online »](#)


Product images:



Immunohistochemistry of paraffin-embedded human gastric cancer using [TA386732] at dilution of 1:100



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue using [TA386732] at dilution of 1:100



Western blot

All lanes: TOB2 antibody at 1.89µg/ml

Lane 1: Mouse skeletal muscle tissue

Lane 2: Mouse heart tissue

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 37, 28 kDa

Observed band size: 37 kDa