

## Product datasheet for **TA386740M**

### **RBMS2 Rabbit Polyclonal Antibody**


#### **Product data:**

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant Human RNA-binding motif, single-stranded-interacting protein 2 protein (1-407AA)
Formulation:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Concentration:	lot specific
Purification:	Antigen Affinity Purified
Conjugation:	Unconjugated
Storage:	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Stability:	1 year from dispatch.
Database Link:	<u><a href="#">Q15434</a></u>




[View online »](#)


## Product images:



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using [TA386740] at dilution of 1:100



Immunohistochemistry of paraffin-embedded human cervical cancer using [TA386740] at dilution of 1:100



Western blot

All lanes: RBMS2 antibody at 0.5µg/ml + MCF-7 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 44 kDa

Observed band size: 44 kDa