

Product datasheet for **TA386773M**

DNAJB8 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200
Reactivity:	Mouse, Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant Human Dnaj homolog subfamily B member 8 protein (1-232AA)
Formulation:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Concentration:	lot specific
Purification:	Antigen Affinity Purified
Conjugation:	Unconjugated
Storage:	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Stability:	1 year from dispatch.
Database Link:	Q8NHS0
Background:	Efficient suppressor of aggregation and toxicity of disease-associated polyglutamine proteins.




[View online »](#)


Product images:



Immunohistochemistry of paraffin-embedded human placenta tissue using [TA386773] at dilution of 1:100



Immunohistochemistry of paraffin-embedded human colon tissue using [TA386773] at dilution of 1:100



Western blot

All lanes: DNAJB8 antibody at 8 μ g/ml

Lane 1: Hela whole cell lysate

Lane 2: Mouse brain tissue

Lane 3: Mouse liver tissue

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 26 kDa

Observed band size: 26 kDa