

Product datasheet for **TA386803M**

HAUS8 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	IF, IHC, WB
Recommended Dilution:	Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200, IF:1:50-1:200
Reactivity:	Mouse, Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant Human HAUS augmin-like complex subunit 8 protein (166-410AA)
Formulation:	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Concentration:	lot specific
Purification:	Antigen Affinity Purified
Conjugation:	Unconjugated
Storage:	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Stability:	1 year from dispatch.
Database Link:	Q9BT25
Background:	Contributes to mitotic spindle assembly, maintenance of centrosome integrity and completion of cytokinesis as part of the HAUS augmin-like complex.




[View online »](#)


Product images:




Immunohistochemistry of paraffin-embedded human kidney tissue using [TA386803] at dilution of 1:100



Immunofluorescent analysis of PC-3 cells using [TA386803] at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)



Immunohistochemistry of paraffin-embedded human tonsil tissue using [TA386803] at dilution of 1:100



Western blot

All lanes: HAUS8 antibody at 2.35 $\mu\text{g/ml}$

Lane 1: Mouse kidney tissue

Lane 2: Mouse brain tissue

Lane 3: A549 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 45, 39 kDa

Observed band size: 45, 50, 28 kDa