

Product datasheet for **TA386807M**

NUP58 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200
Reactivity:	Mouse, Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant Human Nucleoporin p58/p45 protein (211-480AA)
Formulation:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Concentration:	lot specific
Purification:	>95%, Protein G purified
Conjugation:	Unconjugated
Storage:	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Stability:	1 year from dispatch.
Database Link:	Q9BVL2
Background:	Component of the nuclear pore complex, a complex required for the trafficking across the nuclear membrane.




[View online »](#)

Product images:



Immunohistochemistry of paraffin-embedded human small intestine tissue using [TA386807] at dilution of 1:100



Western blot

All lanes: NUP58 antibody at 4 μ g/ml

Lane 1: HepG2 whole cell lysate

Lane 2: HeLa whole cell lysate

Lane 3: LO2 whole cell lysate

Lane 4: Mouse thymus tissue

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 61, 51, 60 kDa

Observed band size: 61 kDa